

## Technical specification

UV penetrability:	1%~100% (1% or 5% steps)
PD:	Single eye: PD:20mm ~ 39.5mm; Double eye: PD: 40mm ~ 85mm (steps 0.05mm)
Lens diameter:	Φ15mm ~ Φ100mm
<b>Display:</b>	<b>5.6" color LED screen</b>
Dot method:	Ink pad
Printer:	High speed thermal printer
Digital connector	RS-232C
Power supply:	AC (10V~240V); 50/60HZ
Max power:	40W
Weight:	3.6kgs
Dimensions:	194mm×210mm×408mm

